
Product Datasheet

Sheep IgM (mu) (orb750733)

Description

Sheep IgM (mu)

Species/Host	Rabbit
Reactivity	Sheep
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Sheep IgM mu heavy chain
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Lyophilized
Concentration	90 mg/mL
Storage	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note	For research use only
Application notes	Anti-Sheep IgM antibody is suitable for use in ELISA, immunohistochemistry, and western blot. Specific conditions for reactivity should be optimized by the end user.
Isotype	Antiserum
Clonality	Polyclonal
Purity	This product was prepared from monospecific antiserum by a delipidation and defibrination. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-rabbit serum, Sheep IgM and Sheep Serum. No reaction was observed against Sheep IgG.
Dilution Range	ELISA: 1:20,000 - 1:100,000, IHC: 1:1,000 - 1:5,000, WB: 1:2,000 - 1:10,000
Expiration Date	12 months from date of receipt.