



Product Datasheet

Rkhd1 Antibody (orb75075)

Description	Rkhd1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q86XN8
Expiration Date	12 months from date of receipt.



Western blot analysis of MDA-MB-361 cell...



Immunohistochemical analysis of formalin...