



Product Datasheet

NIPSNAP3A Antibody (orb75111)

Description	NIPSNAP3A Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9UFN0
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse brain tis...



Immunohistochemical analysis of formalin...