

Product Datasheet

ZBTB7B antibody (orb751123)

Description

Mouse monoclonal antibody to ZBTB7B

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Recombinant full-length human ZBTB7B protein
Target	ZBTB7B
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	200 µg/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG2b
Clonality	Monoclonal
MW	58/80kDa
Uniprot ID	O15156
Entrez	51043
Dilution Range	Flow Cytometry (1-2ug/million cells); Immunofluorescence (1-2ug/ml); Western Blot (1-2ug/ml); ,Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.