
Product Datasheet

ZBTB7B antibody (orb751124)

Description

Mouse monoclonal antibody to ZBTB7B

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, WB

Applications
Immunogen

Recombinant full-length human ZBTB7B protein

Target

ZBTB7B

Preservatives

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

200 µg/ml

Storage

 Antibody with azide - store at 2 to 8 °C.
 Antibody without azide - store at -20 to -80 °C.
 Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Mouse / IgG2b

Clonality

Monoclonal

MW

58/80kDa

Uniprot ID
[O15156](#)
Entrez
[51043](#)
Dilution Range

 Flow Cytometry (1-2ug/million cells);
 Immunofluorescence (1-2ug/ml); Western Blot (1-2ug/ml); ,Optimal dilution for a specific application should be determined.

Expiration Date

12 months from date of receipt.