

---

## Product Datasheet

### PARVG antibody (orb751404)

Evolution Biocomponents.  
**Description**

Rabbit monoclonal antibody to PARVG

<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P
<b>Immunogen</b>	GST-tagged recombinant protein corresponding to human Gamma-parvin
<b>Target</b>	PARVG
<b>Preservatives</b>	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
<b>Concentration</b>	200 µg/ml
<b>Storage</b>	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit / IgG
<b>Clonality</b>	Monoclonal
<b>MW</b>	37kDa
<b>Uniprot ID</b>	<a href="#">Q9HB10</a>
<b>Entrez</b>	<a href="#">64098</a>
<b>Dilution Range</b>	Immunohistochemistry (Formalin-fixed) (1-2ug/ml for 30 minutes at RT),(Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95 &degC followed by cooling at RT for 20 minutes),Optimal dilution for a specific application should be determined.
<b>Expiration Date</b>	12 months from date of receipt.