

Product Datasheet

ZO-1/TJP1 Antibody (orb751451)

Description	ZO-1/TJP1 Antibody
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant C-terminus ZO-1
Form/Appearance	Affinity Purified
Storage	Store at 4°C . Avoid repeated freeze-thaw cycles
Note	For research use only
Isotype	Goat IgG
Uniprot ID	Q07157
Expiration Date	12 months from date of receipt.