

Product Datasheet

TAG-72 (Tumor-Associated Glycoprotein 72) Antibody (orb751514)

Description	TAG-72 (Tumor-Associated Glycoprotein 72) Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant Human TAG-72
Form/Appearance	Affinity Purified
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG
Uniprot ID	Q9XVS1
Entrez	182875
Dilution Range	Immunohistochemistry: 1:1,000-3,000 Western Blot: 1:500-1,000
Expiration Date	12 months from date of receipt.