



Product Datasheet

PCDH12 Antibody (orb75160)

Description	PCDH12 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9NPG4
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 cell lysat...



Immunohistochemical analysis of formalin...