



Product Datasheet

LRFN5 Antibody (orb75164)

Description	LRFN5 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q96N16
Expiration Date	12 months from date of receipt.



Western blot analysis of EL4 cell lysate...



Immunohistochemical analysis of formalin...