

Product Datasheet

5HT (Serotonin) 1E Receptor, 5HT1E R Antibody (orb751641)

Description	5HT (Serotonin) 1E Receptor, 5HT1E R Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic peptide - KLH conjugated 3rd cytoplasmic domain of human 5-HT1E.
Form/Appearance	Affinity Purified
Storage	Aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	P28566
Entrez	3354
Dilution Range	IHC 1:1000-1:2000
Expiration Date	12 months from date of receipt.