

Product Datasheet

MDM2 (phospho Ser186/S188) antibody (orb751959)

Description	Rabbit polyclonal antibody to MDM2 (phospho Ser186/S188)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human MDM2 around the phosphorylation site of Ser186 and Ser188. AA range:151-200
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	/
Uniprot ID	Q00987
Dilution Range	WB: 1: 500-2000, IF: 1: 200-1000, ELISA: 1: 20000
Expiration Date	12 months from date of receipt.