

Product Datasheet

IRAK1 (phospho-Thr387) antibody (orb751960)

Description	Rabbit polyclonal antibody to IRAK1 (phospho-Thr387)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	The antiserum was produced against synthesized peptide derived from human IRAK1 around the phosphorylation site of Thr387. AA range:353-402
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	/
Uniprot ID	P51617
Dilution Range	WB: 1: 500-2000, IF: 1: 200-1000, ELISA: 1: 5000
Expiration Date	12 months from date of receipt.