

Product Datasheet

Bax (phospho-Thr167) antibody (orb751961)

Description	Rabbit polyclonal antibody to Bax (phospho-Thr167)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	The antiserum was produced against synthesized peptide derived from human Bax around the phosphorylation site of Thr167. AA range:133-182
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	/
Uniprot ID	Q07812
Dilution Range	IF: 1: 200-1000, ELISA: 1: 40000
Expiration Date	12 months from date of receipt.