

---

## Product Datasheet

**TGF beta Receptor II (phospho-Ser225/250)  
antibody (orb751977)**

**Description**

Rabbit polyclonal antibody to TGF

Ser225/250)

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse

**Conjugation**

Unconjugated

**Tested Applications**

ELISA, IF, IHC-P

**Immunogen**

The antiserum was produced against synthesized peptide derived from human TGF beta Receptor II around the phosphorylation site of Ser225/250. AA range:191-240

**Preservatives**

Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.

**Concentration**

1mg/ml

**Storage**

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note**

For research use only

**Clonality**

Polyclonal

**MW**

~ 65 kDa

**Uniprot ID****P37173****Dilution Range**

WB: 1: 500-2000, ICC: 1: 100-300,  
ELISA: 1: 10000

**Expiration Date**

12 months from date of receipt.