

Product Datasheet

PLK1 (phospho-Thr210) antibody (orb751978)

Description	Rabbit polyclonal antibody to PLK1 (phospho-Thr210)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human PLK1 around the phosphorylation site of Thr210. AA range:176-225
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 70 kDa
Uniprot ID	P53350
Dilution Range	WB: 1: 500-2000, ICC: 1: 100-300, ELISA: 1: 5000
Expiration Date	12 months from date of receipt.