

## Product Datasheet

### IGFBP-3 (phospho-Ser183) antibody (orb751996)

<b>Description</b>	Rabbit polyclonal antibody to IGFBP-3 (phospho-Ser183)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, IHC-P, WB
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from human IGFBP-3 around the phosphorylation site of Ser183. AA range:151-200
<b>Preservatives</b>	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>MW</b>	~ 31 kDa
<b>Uniprot ID</b>	<b>P17936</b>
<b>Dilution Range</b>	WB: 1: 500-2000, ICC: 1: 100-300, ELISA: 1: 5000
<b>Expiration Date</b>	12 months from date of receipt.