

Product Datasheet

AurB/C (phospho-Thr236/202) antibody (orb752002)

Description	Rabbit polyclonal antibody to AurB/C (phospho-Thr236/202)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P
Immunogen	The antiserum was produced against synthesized peptide derived from human AurB/C around the phosphorylation site of Thr236/202. AA range:201-250
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	/
Uniprot ID	Q9UQB9
Dilution Range	ICC: 1: 100-300, ELISA: 1: 10000
Expiration Date	12 months from date of receipt.