

Product Datasheet

ITIH4(70k,Cleaved-Arg661) antibody (orb752023)

Description	Rabbit polyclonal antibody to ITIH4(70k,Cleaved-Arg661)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P
Immunogen	The antiserum was produced against synthesized peptide derived from human ITIH4. AA range:612-661
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	/
Uniprot ID	Q14624
Dilution Range	ICC: 1: 100-300, ELISA: 1: 10000
Expiration Date	12 months from date of receipt.