

Product Datasheet

PAK4/5/6 (phospho-Ser474) antibody (orb752027)

Description	Rabbit polyclonal antibody to PAK4/5/6 (phospho-Ser474)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human PAK4/5/6 around the phosphorylation site of Ser474. AA range:441-490
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	/
Uniprot ID	O96013
Dilution Range	WB: 1: 500-2000, ICC: 1: 100-300, ELISA: 1: 20000
Expiration Date	12 months from date of receipt.