

Product Datasheet

Lamin A (phospho-Ser22) antibody (orb752086)

Description	Rabbit polyclonal antibody to Lamin A (phospho-Ser22)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Lamin A/C around the phosphorylation site of Ser392. AA range:361-410
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 74 kDa
Uniprot ID	P02545
Dilution Range	WB: 1: 500-2000, ICC: 1: 100-300, IF: 1: 200-1000, ELISA: 1: 10000
Expiration Date	12 months from date of receipt.