

Product Datasheet

Claudin 7 (phospho-Tyr210) antibody (orb752151)

Description	Rabbit polyclonal antibody to Claudin 7 (phospho-Tyr210)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Claudin 7 around the phosphorylation site of Tyr210. AA range:162-211
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 32 kDa
Uniprot ID	O95471
Dilution Range	WB: 1: 500-2000, ELISA: 1: 5000
Expiration Date	12 months from date of receipt.