



Product Datasheet

Seasonal H1N1 Hemagglutinin Antibody (orb75223)

Description	Seasonal H1N1 Hemagglutinin Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	BOVX44
Expiration Date	12 months from date of receipt.



Western blot analysis of recombinant (la...



Bar chart illustrates about the ELISA re...