

## Product Datasheet

# CD227/MUC1 (phospho-Tyr1229) antibody (orb752239)

<b>Description</b>	Rabbit polyclonal antibody to CD227/MUC1 (phospho-Tyr1229)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, IHC-P, WB
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from human CD227/MUC1 around the phosphorylation site of Tyr1229. AA range:1201-1250
<b>Preservatives</b>	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>MW</b>	~ 170 kDa
<b>Uniprot ID</b>	<a href="#">P15941</a>
<b>Dilution Range</b>	WB: 1: 500-2000, ICC: 1: 100-300, IF: 1: 200-1000, ELISA: 1: 10000
<b>Expiration Date</b>	12 months from date of receipt.