



## Product Datasheet

# Swine H1N1 Hemagglutinin Antibody (orb75224)

<b>Description</b>	Swine H1N1 Hemagglutinin Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Storage</b>	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<b>C4AL34</b>
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of recombinant (la...



Bar chart illustrates about the ELISA re...