



Product Datasheet

Seasonal H1N1 Hemagglutinin Antibody (orb75227)

| | |
|----------------------------|---|
| Description | Seasonal H1N1 Hemagglutinin Antibody |
| Species/Host | Rabbit |
| Reactivity | Virus |
| Conjugation | Unconjugated |
| Tested Applications | ELISA |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | B0VX44 |
| Expiration Date | 12 months from date of receipt. |



Bar chart illustrates about the ELISA re...



Bar chart illustrates about the ELISA re...