

Product Datasheet

Canine PD-1 / PDCD1 Protein (orb752324)

Description	Canine PD-1, His Tag is expressed from human 293 cells (HEK293). It contains AA Leu 25 - Gly 168 (Ac...
Reactivity	Canine
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PD-1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 18.0 kDa. The protein migrates as 30-40 kDa under reducing (R) condition (SDS-PAGE) due to Glycosylation.
Protein Sequence	NP_001301026.1
Purity	90%
Source	Canine PD-1, His Tag (orb752324) is expressed from human 293 cells (HEK293). It contains AA Leu 25 - Gly 168 (Accession # NP_001301026.1).
MW	18.0 kDa
Uniprot ID	A0A090BAM7
NCBI	NP_001301026.1
Expiration Date	6 months from date of receipt.