


Product Datasheet

Human IL-1 beta / IL-1F2 Protein (orb752327)

Description	Human IL-1 beta, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Ala 117 - Ser 269...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-1 beta
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	N-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the N-terminus. The protein has a calculated MW of 43.8 kDa. The protein migrates as 50-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human IL-1 beta, Fc Tag (orb752327) is expressed from human 293 cells (HEK293). It contains AA Ala 117 - Ser 269 (Accession # P01584-1).
MW	43.8 kDa
Uniprot ID	P01584
NCBI	NP_000567.1
Expiration Date	6 months from date of receipt.



Immobilized
Human IL-1
RII, His Tag
(Cat...