
Product Datasheet

**Mouse / Rat TGF-Beta 1 / TGFB1 Protein
(orb752336)**

Description

Mouse / Rat TGFB1, Tag Free is

(HEK293). It contains AA Ala 279 - Ser...

Reactivity	Mouse, Rat
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TGF-beta 1
Preservatives	50 mM HAC, pH3.0
Form/Appearance	Powder
Storage	-20°C
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 12.8 kDa. The protein migrates as 13 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	NP_035707.1
Purity	95%
Source	Mouse / Rat TGFB1, Tag Free (orb752336) is expressed from human 293 cells (HEK293). It contains AA Ala 279 - Ser 390 (Accession # P04202-1).
MW	12.8 kDa
Uniprot ID	P04202
NCBI	NP_035707.1
Expiration Date	6 months from date of receipt.