


Product Datasheet

Seasonal H1N1 Nonstructural Protein 1 Antibody (orb75258)

Description	Seasonal H1N1 Nonstructural Protein 1 Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	A4U724
Expiration Date	12 months from date of receipt.



Bar chart illustrates about the ELISA re...