

Product Datasheet

KCNJ12 antibody (FITC) (orb752594)

Description	KCNJ12 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human ATP-sensitive inward rectifier potassium channel 12 protein (182-433aa)
Target	KCNJ12
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q14500
Expiration Date	12 months from date of receipt.