

Product Datasheet

OR9G1 antibody (orb752660)

Description	OR9G1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Olfactory receptor 9G1 protein (290-304aa)
Target	OR9G1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NH87
Dilution Range	IF: 1:100-1:500
Expiration Date	12 months from date of receipt.