

Product Datasheet

SMARCC2 antibody (FITC) (orb752749)

Description	SMARCC2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human SWI/SNF complex subunit SMARCC2 protein (300-647aa)
Target	SMARCC2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8TAQ2
Expiration Date	12 months from date of receipt.