

Product Datasheet

Mouse mEPOA Protein (orb752934)

Description	EPO Mouse Recombinant produced in Baculovirus is a single glycosylated polypeptide chain containing ...
Preservatives	EPO Mouse protein (0.25mg/ml) contains 10% glycerol and Phosphate-Buffered Saline (pH 7.4).
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Cytokines And Growth Factors
Protein Sequence	ADPMAPPRLI CDSRVLERYI LEAKEAENV MGCAEGPRLS ENITVPDTKV NFYAWKRMEV EEQAIEVWQG LSELLSEAILQ AQALLANSSQ PPETLQLHID KAISGLRSLT SLLRVLGAQK ELMSPDTP PAPLRTLTVDFCKLFRVYA NFLRGKLLY TGEVCRRGDR HHHHHH
Purity	Greater than 95.0% as determined by SDS-PAGE.
Source	Sf9, Baculovirus cells
Expiration Date	12 months from date of receipt.