

Product Datasheet

Porcine pGMCSF Protein (orb752942)

Description	GMCSF Poricne Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain contain...
Preservatives	GMCSF protein solution (0.5mg/ml) containing Phosphate-Buffered Saline (pH7.4).
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Cytokines And Growth Factors
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MAPTRPPSPV TRPWQHVD AI KEALSLNNS NDTAAVMNET VDVVCEMFDP QEPTCVQTRL NLYKQGLRGS LTRLKSPLTL LAKHYEQHCP LTEETSCETQ SITFKSFKDS LNKFLFTIPF DCWGPVKK
Purity	Greater than 95.0% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.