
Product Datasheet

Human ALPP Protein (orb753009)

Description

ALPP Human produced in Sf9 Baculovirus cells

containing ...

Preservatives

ALPP protein solution (0.5mg/ml) containing Phosphate-Buffered Saline (pH 7.4) and 10% glycerol.

Form/Appearance

Sterile Filtered colorless solution.

Storage

Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles

Note

For research use only

Application notes

Enzymes

Protein Sequence

ADLMIIPVEE ENPDFWNREA AEALGAACKL
 QPAQTAANKL IIFLGDGMGV STVTAARILK
 GQKKDKLGPE LPLAMDRFPY VALSKTYNVD
 KHVPDSGATA TAYLCGVKGN FQTIGLSAAA
 RFNQCNTTRG NEVISVMNRA KKAGKSVGUV
 TTTRVQHASP AGTYAHTVNR NWYSDADVPA
 SARQEGCQDI ATQLISNMDI DVLGGGRKY
 MFRMGTPDPE YPDDYSQGGT RLDGKNLVQE
 WLAKRQGARY VWNRTLMQA SLDPSVTHLM
 GLFEPGDMKY EIHRDSTLDP SLMEMTEAAL
 RLLSRNPRGF FLFVEGGRID HGHHESTRAYR
 ALTETIMFDD AIERAGQLTS EEDTSLVTA
 DSHVFSFGG YPLRGSSIFG LAPGKARDRK
 AYTVELLYGNG PGYVLKDGAR PDVTESESGS
 PEYRQQSAVP LDEETHAGED VAVFARGPQA
 HLVHGVQEQT FIAHVMAFAA CLEPYTACDL
 APPAGTTDHH HHHH

Purity

Greater than 95.0% as determined by SDS-PAGE.

Source

Sf9, Baculovirus cells

Expiration Date

12 months from date of receipt.