

Product Datasheet

Human BACE1 Protein (orb753010)

Description	BACE1 Human produced in Sf9 Baculovirus cells is a single, glycosylated polypeptide chain containing...
Preservatives	BACE1 protein solution containing Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Enzymes
Protein Sequence	TQHGIRLPLR SGLGGAPLGL RLPRETDEEP EEPGRRSFV EMVDNLRGKS GQGYVEMTV GSPPQTLNIL VDTGSSNFAV GAAPHPFLHR YYQRQLSSTY RDLRKGVVYP YTQGWEGEL GTDLVSIPIHG PNVTVRANIA AITESDKFFI NGSNWEGILG LAYAEIARPD DSLEPPFDSL VKQTHVPLNF SLQLCGAGFP LNQSEVLASV GGSMIIGGID HSLYTGSLWY TPIRREWYEE VIIVRVEING QDLKMDCKEY NYDKSIVDSG TTNLRPLPKV FEAAVKSIA ASSTEKFPDG FWLGEQLVCW QAGTTPWNIF PVISLYLMGE VTNQSFRTI LPQQYLRPVE DVATSQDDCY KFAISQSSTG TVMGAVIMEG FYVVFDRARK RIGFAVSACH VHDEFRTAAV EGPVFTLDME DCGYNIPQTD ESTLMTHHHH HH
Purity	Greater than 90.0% as determined by SDS-PAGE.
Source	Sf9, Baculovirus cells
Expiration Date	12 months from date of receipt.