



Product Datasheet

OCIAD2 Antibody (orb75303)

| | |
|----------------------------|---|
| Description | OCIAD2 Antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IF, WB |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | Q56VL3 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of SK-N-SH cell ly...



Confocal immunofluorescence analysis of ...