
Product Datasheet

Human CD72 Protein (orb753100)

Description

CD72 Human Recombinant produced in Sf9

polypeptide chain...

Preservatives

CD72 protein solution (0.5mg/ml) contains 50mM Tris-HCl buffer (pH6.8), 0.2M NaCl, 2mM DTT and 50% glycerol.

Form/Appearance

Sterile Filtered colorless solution.

Storage

Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles

Note

For research use only

Application notes

Recombinant & Natural Proteins

Protein Sequence

ADPRYLQVSQ QLQQTNRVLE VTNSSLRQQL
 RLKITQLGQS AEDLQGSRRE
 LAQSQEALQVEQRAHQAAEG QLQACQADRQ
 KTKETLQSEE QRRALEQKL SNMENRLKPF
 FTCGSADTCCPSGWIMHQKS CFYISLTSKN
 WQESQKQCET LSSKLATFSE IYPQSHSYF
 LNSLLPNGGSGNSYWTGLSS NKDWKLTDDT
 QRTRTYAQSS KCNKVHKTWS WWTLESESCR
 SSLPYICEMTAFRFPDLEPK SCDKTHTCPP
 CPAPELLGGP SVFLFPPKPK DTLMISRTPE
 VTCVVVDVSHEDPEVKFNWY VDGVEVHNAK
 TKPREEQYNS TYRVVSVLTV LHQDWLNGKE
 YKCKVSNKALPAPIEKTISK AKGQPREPQV
 YTLPPSRDEL TKNQVSLTCL VKGFYPSDIA
 VEWESNGQPENNYKTTTPVL DSDGSFFLYS
 KLTVDKSRWQ QGNVFSCSVM HEALHNHYTQ
 KSLSLSPGKH HHHHH

Purity

Greater than 90% as determined by SDS-PAGE.

Source

Sf9, Baculovirus cells

Expiration Date

12 months from date of receipt.