

Product Datasheet

Human PLCXD3 Protein (orb753157)

Description	PLCXD3 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containi...
Preservatives	The PLCXD3 solution (0.25mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 20% glycerol, 0.15M NaCl and 1mM DTT.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGSMASSQGK NELKLADWMA TLPESMHSIP LTNLAIPGSH DSFSFYIDEA SPVGPEQPET VQNFVSVFGT VAKKLMRKWL ATQTMNFTGQ LGAGIRYFDL RISTKPRDPD NELYFAHGLF SAKVNEGLEE INAFILDHKK EVVFLDFNHF YGMQKYHHEK LVQMLKDIYG NKMCPAIFAQ EVSLKYLWEK DYQVLVIFYHS PVALEVPFLW PGQMMPAPWA NTTDPEKLIQ FLQASITERR KKGSSFFISQV VLTPKASTVV KGVASGLRET ITERALPAMM QWVRTQKPGE SGINIVTADF VELGDFISTV IKLNYVFDEG EANT
Purity	Greater than 85.0% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.