

Product Datasheet

Human SELPLG Protein (orb753159)

Description	Recombinant Human Selectin P Ligand, HEK
Preservatives	The SELPLG solution (1mg/1ml) contains phosphate buffered saline (pH7.4) and 10% glycerol.
Form/Appearance	Sterile filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	DGSQATEY EY LDYDFLPETE PPEMLRNSTD TTPLTGPGTP ESTTVEPAAR RSTGLDAGGA VTELTTELAN MGNLSTDSAA MEIQTTQPAA TEAQTTPPLAA TEAQTTRLTA TEAQTTPPLAA TEAQTTPPAA TEAQTTOPTG LEAQTTPAA MEAQTTPAA MEAQTTPPAA MEAQTQTTA MEAQTTAPEA TEAQTTOPTA TEAQTTPPLAA MEALSTEPSA TEALSMEPTT KRGLFIPFSV SSVTHKGIPM AASNLSVLEP KSCDKTHTCP PCPAPELLGG PSVFLFPPKP KDTLMISRTP EVTCVVVDVS HEDPEVKFNW YVDGVEVHNA KTKPREEQYN STYRVVSVLT VLNQDNLNGK EYKCKVSNKA LPAPIEKTIS KAKGQPREPQ VYTLPPSRDE LTKNQVSLT LVKGFYPSDI AVEWESNGQP ENNYKTTTPV LDSDGSFFLY SKLTVDKSRW QQGNVFCVSV MHEALHNHYT QKLSLSLSPGK HHHHHH
Purity	Greater than 90.0% as determined by SDS-PAGE.
Source	HEK293 Cells
Expiration Date	12 months from date of receipt.