

---

## Product Datasheet

### Human SIGLEC10 Protein (orb753163)

**Description**

SIGLEC10 Human produced in Sf9 Insect cells

containing 6...

**Preservatives**

The SIGLEC10 solution (0.25mg/ml) contains Phosphate Buffered Saline (pH 7.4) and 10% glycerol.

**Form/Appearance**

Sterile Filtered colorless solution.

**Storage**

Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles

**Note**

For research use only

**Application notes**

Recombinant &amp; Natural Proteins

**Protein Sequence**

```

MDGRFWIRVQ ESVMVPEGLC ISVPCSFSSYP
RQDWTGSTPA YGYWFKAVTE
TTKGAPVATNHQSREVEMST RGRFQLTGDP
AKGNCSLVIR DAQMQDESQY FFRVERGSYV
RYNFMNDGFFLKVTALTQKP DVYIPETLEP
GQPVTVICVF NWAEECPFP SFSWTGAALS
SQGTKPTTSHFVLSFTPRP QDHNTDLTCH
VDFSRKGVSA QRTVRLRVAY APRDLVISI
RDNTPALEPQPQGNVPYLEA QKGQFLRLLC
AADSQPPATL SWVLQNRVLS SSHPWGPRPL
GLELPGVKAGDSGRYTCRAE NRLGSQQRAL
DLSVQYPPEN LRVMVSQANR TVLENLGNGT
SLPVLEGQSLCLVCVTHSSP PARLSWTQRG
QVLSPSQPSD PGVLELPRVQ VEHEGEFTCH
ARHPLGSQHVLSLSVHYKK GLISTAFSNL
EPKSCDKTHT CPPCPAPELL GGPSVFLFPP
KPKDTLMISRTPEVTCVVVD VSHEDPEVKF
NWYVDGVEVH NAKTKPREEQ YNSTYRVVSV
LTVLHQDWLNGKEYKCKVSN KALPAPIEKT
ISKAKGQPRE PQVYTLPPSR DELTKNQVSL
TCLVKGFYPSDIAVEWESNG QPENNYKTPP
PVLDSGDGFF LYSKLTVDKS RWQQGNVFSC
SVMHEALHNNH YTKQSLSLSP GKHHHHHHH
  
```

**Purity**

Greater than 85.0% as determined by SDS-PAGE.

**Source**

Sf9, Baculovirus cells

**Expiration Date**

12 months from date of receipt.