

Product Datasheet

Mouse mSOST Protein (orb753168)

Description	SOST Mouse Recombinant produced in HEK293 cells is a single, glycosylated polypeptide chain (24-211 ...
Preservatives	SOST protein (0.25mg/ml) contains 10% glycerol and Phosphate-Buffered Saline (pH 7.4).
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	QGWQAFRND A TEVIPGLGEY PEPPPENNT MNRAENGGR P P H P Y D A K D V S E Y S C R E L H Y T R F L T D G P C R S A K P V T E L V C S G Q C G P A R L L P N A I G R V K W W R P N G P D F R C I P D R Y R A Q R V Q L L C P G G A A P R S R K V R L V A S C K C K R L T R F H N Q S E L K D F G P E T A R P Q K G R K P R P G A R G A K A N Q A E L E N A Y H H H H H H
Purity	Greater than 90.0% as determined by SDS-PAGE.
Source	HEK293 cells
Expiration Date	12 months from date of receipt.