

Product Datasheet

CoV Spike antibody (orb753252)

Description	Mouse monoclonal antibody to CoV Spike
Reactivity	Virus
Conjugation	Unconjugated
Immunogen	The antibody was developed using purified bovine coronavirus. (eptopoe not mapped)
Preservatives	1mg/ml in 1X PBS pH-7.2 & 0.01% NaN ₃ .
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Monoclonal Antibodies, Anti-Viral
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.