

Product Datasheet

CoV Spike antibody (orb753254)

Description	Mouse monoclonal antibody to CoV Spike
Reactivity	Virus
Conjugation	Unconjugated
Immunogen	The antibody was developed using purified porcine coronavirus clone PTT1E11. (eptopoe not mapped)
Preservatives	1mg/ml in 1X PBS pH-7.2 & 0.01% NaN ₃ .
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Applications: This porcine Coronavirus Spike monoclonal antibody can be used in WB and ELISA
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.