

Product Datasheet

Rabbit CoV-2 S1 antibody (orb753260)

Description	Mouse monoclonal antibody to Rabbit CoV-2 S1
Reactivity	Virus
Conjugation	Unconjugated
Immunogen	The native monoclonal antibody was generated by sequencing peripheral blood lymphocytes of a patient exposed to the SARS-CoV
Preservatives	1mg/ml in PBS and 0.02% Proclin 300.
Form/Appearance	Sterile Filtered liquid formulation.
Storage	Stability: Store at 4°C, stable for 2 weeks. For long-term storage, store at -20°C
Note	For research use only
Application notes	Applications: ELISA
Clonality	Polyclonal
Source	HEK293 Cells
Expiration Date	12 months from date of receipt.