

## Product Datasheet

### SARS Nucleocapsid antibody (Biotin) (orb753262)

<b>Description</b>	Mouse monoclonal antibody against SARS Nucleocapsid conjugated to Biotin
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Immunogen</b>	Recombinant protein fragment 1-49 a.a. of the SARS nucleoprotein
<b>Preservatives</b>	1mg/ml in DPBS, 1% rAlbumin and 0.1% Proclin 950.
<b>Form/Appearance</b>	Sterile Filtered liquid formulation.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Applications: ELISA
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE gel.
<b>Expiration Date</b>	12 months from date of receipt.