

Product Datasheet

SARS Nucleocapsid antibody (PE) (orb753263)

Description	Mouse monoclonal antibody against SARS Nucleocapsid conjugated to PE
Reactivity	Virus
Conjugation	Unconjugated
Immunogen	Recombinant protein fragment 1-49 a.a. of the SARS nucleoprotein
Preservatives	1mg/ml in DPBS, 1% rAlbumin and 0.1% Proclin 950.
Form/Appearance	Sterile Filtered liquid formulation.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Applications: ELISA
Clonality	Monoclonal
Purity	Greater than 90% as determined by SDS-PAGE gel.
Expiration Date	12 months from date of receipt.