

## Product Datasheet

### SARS MERS Spike S1 antibody (orb753264)

<b>Description</b>	Mouse monoclonal antibody to SARS MERS Spike S1
<b>Reactivity</b>	Virus
<b>Preservatives</b>	1mg/ml in PBS and 0.1% NaN <sub>3</sub> .
<b>Form/Appearance</b>	Sterile Filtered liquid formulation.
<b>Storage</b>	Stability: Store at 4°C, stable for 2 weeks. For long-term storage, store at -20°C
<b>Note</b>	For research use only
<b>Application notes</b>	Applications: ELISA
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE gel.
<b>Expiration Date</b>	12 months from date of receipt.