

Product Datasheet

SARS MERS Spike S2 antibody (orb753265)

Description	Mouse monoclonal antibody to SARS MERS Spike S2
Reactivity	Virus
Conjugation	Unconjugated
Immunogen	Recombinant MERS-CoV Spike S2 Subunit (752-1296aa) purified from Baculovirus
Preservatives	1mg/ml containing Phosphate-Buffered Saline (pH 7.4), 0.02% Sodium Azide and 10% Glycerol.
Form/Appearance	Sterile filtered colorless solution
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Applications: SARS MERS Spike S2 antibody has been tested by ELISA and Western blot analysis to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results
Isotype	IgG2b kappa
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.